# NISSEI BLOOD PRESSURE MONITOR DSK-1031

## INSTRUCTIONS

This manual is intended to assist you in the safe and efficient operation of BLOOD PRESSURE MONITOR DSK-1031. The product must be used in accordance with the procedures contained in this manual and must not be used for purposes other than those described herein. It is important to read and understand the entire manual. In particular, please read carefully and become familiar with the section entitled "TIPS ON TAKING YOUR BLOOD PRESSURE".

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
PL

# CONTENTS

<b>PART NAMES AND PRODUCT COMPONENTS</b> .....	2	<b>MAKING A MEASUREMENT</b> .....	8
<b>GENERAL INFORMATION</b> .....	2	<b>VIEWING THE SAVED READINGS (MEMORY FUNCTION)</b> ..	11
<b>PRECAUTIONS AND CONTRAINDICATIONS</b> .....	5	<b>TROUBLESHOOTING</b> .....	12
<b>TIPS ON TAKING YOUR BLOOD PRESSURE</b> .....	6	<b>TECHNICAL DESCRIPTION FOR ELECTROMAGNETIC DISTURBANCES</b> .....	13
<b>PREPARATION BEFORE USE</b> .....	7		

## PART NAMES AND PRODUCT COMPONENTS

## GENERAL INFORMATION

### INDICATIONS FOR USE

This product is intended for noninvasive measurement of systolic and diastolic blood pressure, determination of pulse rate and calculation of pulse pressure in adults in a home healthcare environment. The product is not designed for neonatal use. Please consult with


your doctor or physician to use this product to take blood pressure of child or person in pregnancy or under pre-eclamptic condition.

### METHOD OF MEASUREMENT

This product employs the oscillometric method for measurement of blood pressure and pulse rate. The cuff is connected to the main unit and wrapped around the arm. Circuits within the cuff sense the small oscillations in pressure against the cuff produced by the expansion and contraction of the arteries in the arm in response to each heart beat. The amplitude of each pressure waves is measured, converted to millimeters of mercury, and displayed on the LCD as a digital value.

### FUZZY INFLATION

This product employs FUZZY inflation to ease stress on arm during cuff inflation. By detecting pulse signals while inflating, this product is able to determine the sufficient pressure for each measurement.

### CUFF SYMBOL

This product displays whether the cuff is applied properly or not by analyzing the pressure wave. [ (OK) ] indicates that the cuff is applied properly but [ (O) ] indicates that the cuff may be wrapped too tightly or too loosely.

### RELIABILITY SYMBOL

Some of the important aspects in measuring blood pressure are how properly the cuff is applied and how calm you are during measurement. These aspects can affect measurement values and such values cannot said to be correct values. This product shows [ (G) ] when the measurement is supposedly made in good condition without influences from such aspects.

### WHO CLASSIFICATION SYMBOL

Measured blood pressure values are classified against WHO guideline (1999). The higher class applies when systolic and diastolic pressures fall into different classes. For specific information on your blood pressure, contact your physician. Never make any judgment on your own based on the classification indication only.

DISPLAY	WHO CLASSIFICATION	SYS	DIA
	Hypertension (severe)	>180	>110
	Hypertension (moderate)	160-179	100-109
	Hypertension (mild)	140-159	90-99
	High Normal	130-139	85-89
	Normal	120-129	80-84
	Optimal	<120	<80

### PULSE PRESSURE DISPLAY

This product calculates and displays pulse pressure.


While systolic blood pressure keeps increasing with age, diastolic blood pressure tends to start decreasing at around age 50. Pulse pressure is the value shown when the diastolic is subtracted from the systolic measurement and therefore pulse pressure tends to increase with age. Pulse pressure has been considered to be related to stiffness of artery and has been studied as one of circulatory risk factors. It is sometimes said that 45 mmHg of pulse pressure is normal value. Arterial sclerosis cannot be represented by pulse pressure only. However, observation of pulse pressure in long term will be important.


### BODY MOTION SYMBOL

Blood pressure value taken while moving cannot be said to be the correct value because body movement can affect blood pressure.

This product analyzes pulse wave and displays [ (D) ] when body motion is detected. [ (D) ] indicates the results might be affected by body movement.

## IRREGULAR PULSE RHYTHM SYMBOL

Pulse rhythm can be disturbed from talking, moving or arrhythmias. This product displays [  ] when the difference of shortest and longest time intervals between pulse waves is over 25%, indicating irregular pulse rhythm.

An example of irregular pulse rhythm compared to regular pulse rhythm

## MEMORY FUNCTION

The measured values are automatically saved for later review in either one of two memory banks, [1] or [2]. These two memory banks can be used to save readings of two persons separately or to save morning readings and evening readings separately.

Each bank can save up to 60 readings. When the number of stored readings reaches 60, the oldest reading will be deleted to record new reading.

Saved readings are recalled with measurement date and time when clock function is activated. \*The clock needs to be activated before measurement for date and time to be saved with the measured valued.

[Err] results are not saved.

## CARE AND MAINTENANCE

Keep the product clean. Inspect its cleanliness after use. To clean, use only a soft dry cloth. Do not use gasoline, paint thinner, or other strong solvents. Since the cuff may absorb perspiration and other fluids, inspect it for stain and discoloration after each use. When cleaning the cuff, use a synthetic detergent and gently rub the surface. Air dry thoroughly. Make sure fluid never gets in AIR HOSE. Do not machine wash or scrub it.

When storing the product, do not place heavy objects on it. Do not fold or bend the cuff forcibly. To disconnect the cuff from the main unit, do not pull AIR PLUG but hold AIR PLUG and unplug it. Do not coil AIR HOSE too tightly. In case the product is stored in the environment with ambient temperature above 40°C or below 10°C, please leave it for at least 2 hours before taking a measurement. Remove the batteries if the product is to be stored for an extended period of time. Keep the batteries out of reach of children.

Stop using the product immediately and contact your dealer or the manufacturer in case any

visible damage is found on the device.

We suggest that you have your blood pressure monitor checked every 2 years. This operation may only be performed by the manufacturer or by firms authorized by the manufacturer.

## PRODUCT SPECIFICATIONS

Model	: DSK-1031
Operating principle	: Oscillometric method
Indicator	: 15 digits liquid crystal display
Pressure indicating range	: 0 to 300 mmHg (cuff pressure)
Measuring range	: 50 to 250 mmHg (systolic), 40 to 180 mmHg (diastolic), 40 to 160 bpm (pulse rate)
Accuracy*	: $\pm 3$ mmHg (cuff pressure), $\pm 5\%$ of reading (pulse rate)
Inflation	: Automatic inflation (FUZZY)
Deflation	: Automatic (electric control valve)
Exhaust	: Automatic exhaust valve
Power supply	: Four 1.5 volt LR6 (AA alkaline) batteries or AC adaptor ADP-W5 series
Electric rating	: DC6V/4W, with four LR6 batteries Adaptor: AC100-240V, 50-60Hz, 0.12A, Unit: DC6V, $\equiv$ 500mA, with designated AC adaptor
Memory	: 2 banks, each stores 60 readings
Operating condition	: +10°C to +40°C, 15% to 85% RH (noncondensing)
Transportation/storage condition	: -20°C to +60°C, 10% to 95% RH (noncondensing)
Cuff	: Model; DSC-EP10, Coverage arm circumference; 22 to 42 cm, applied part; BF
Main unit	: Weight; Approx. 250 g (without batteries) Size; Approx. 115 x 115 x 65.9 mm (W x D x H)
Protection class IP	: IP 20: Protected against solid foreign particles with a diameter of more than 12.5 mm, no protection against water.
Protection against electric shock	: Internally powered equipment/Class II equipment, Type BF applied part

Mode of operation	: Continuous operation
Classification	: Class II / Internally powered equipment
Expected service life	: 5 years
Key to symbols	 : Type BF applied part  : Class II equipment  : Refer to instruction manual/booklet.  : Keep dry  : The used electrical and electronic products are not household waste. Follow your national/local recycling rules to dispose of them properly. In the EU countries, please refer to waste management symbol(s) marked on the package or the instrument.

This device complies with EN1060-1:1995+A2:2009 Non-invasive sphygmomanometers Part 1: General requirements and EN1060-3:1997+A2:2009 Non-invasive sphygmomanometers Part 3: Supplementary requirements for electro-mechanical blood pressure measuring system.

\*Accuracy is guaranteed with the measured values that are within the measuring range.

Specifications are subject to change without notice due to improvements in performance.

\*The measurement accuracy of the product has been proven according to ISO 81060-2 protocol. In the clinical study, K5 was used for the determination of diastolic pressure values at all auscultatory measurements.

\*This product is intended for use in the environment with one atmospheric pressure.

## PRECAUTIONS AND CONTRAINDICATIONS

Do not use this product without consultation with your doctor if you are under dialysis therapy or on anticoagulants, antiplatelets or steroids. Use of this instrument under such conditions could cause internal bleeding.

For specific information on your blood pressure, contact your physician. Never make any judgment on your own regarding measurement results.

Do not use this product along with implantable and wearable medical electrical equipment such as pacemaker, defibrillator, or electrocardiographic monitor. This product is not also intended to be used with HF surgical equipment.

Do not use this product in an explosive environment such as near flammable anesthetics or inside oxygen chamber.

The system may fail to yield specified measurement accuracy if operated or stored in temperature or humidity conditions outside the limits stated in the specifications section of this manual.

Do not use cuffs or accessories other than those specified by the manufacturer.

Otherwise, correct measurement readings cannot be obtained.

Use of this product in areas near mobile phones, microwave ovens or other devices with strong electromagnetic field may cause malfunctions.

Do not apply the cuff over wounded arm, arm under an intravascular access or therapy or an arterio-venous shunt, or arm on the side of a mastectomy or lymph node clearance. Otherwise injury may be resulted.

Make sure that inflation of the cuff is not causing prolonged impairment of blood circulation. Also, be cautious about temporary loss of the functions of any other medical equipment if any monitoring equipment is used on the same limb with the blood pressure measuring cuff.

To avoid harmful injury due to interfered blood flow from cuff inflation,

- Make sure that AIR HOSE is not kinking before measurement. Otherwise, cuff inflation may not be conducted properly and prolonged, and
- Do not make measurements repeatedly.

To avoid any possibility of accidental strangulation, keep this product away from children and do not drape AIR HOSE around your neck.

Because the product includes precision parts, avoid extreme temperature variations, humidity, shock, dust, lint, and direct sunlight. Do not drop or strike the product. Make sure not to expose it to moisture. This product is not water resistant.

Do not press the display and operation keys or place the monitor with display face down.

Do not take out batteries or unplug the AC adaptor when the monitor is turned on. Make sure to switch off the monitor before removing batteries or AC adaptor.

Do not touch the output plug of AC adaptor during measurement.

Do not disassemble or modify the product.

## TIPS ON TAKING YOUR BLOOD PRESSURE

Blood pressure is the force exerted by the heart in pumping the blood through the arteries and the resistance by the veins in this flow.

**Blood pressure varies all the time, influenced by mental and physical factors, and is never constant.**

In general, blood pressure is highest during the working hours and gradually decreases during the afternoon and evening hours. It is low during sleep and increases at a relatively fast rate after arising from bed.

Causes for Changes in Blood Pressure

- Body movement
- Emotions
- Recent Urination or Bowel Movement
- Changes in the surroundings such as movement or noise, etc.
- Conversation
- Eating
- Mental Tension
- Drinking Alcohol
- Room temperature
- Nervousness
- Smoking

**Blood pressure measured at home tends to be lower than when measured in a hospital, clinic or doctor's office.**

Do not inflate the cuff when it is not wrapped around your arm.

The device contains small parts and batteries which could be swallowed by children or pets. They should therefore be kept out of the reach of children and pets at all times.

This product is not designed for self-use by unspecified persons in public areas.

Do not apply the cuff on the limb which the intravenous drip infusion is implemented.

Any serious incident occurred in relation to the device should be reported to the manufacturer and the competent authority in your country/area. If you have no contact information of such authority, please contact the manufacturer or EU authorized representative whose contact information is indicated in this instruction manual.

This is because you are tense at the hospital and relaxed at home. It is important to know your stable normal blood pressure at home.

**Let a qualified physician interpret your blood pressure readings.**

Depending on your age, weight and general condition, blood pressure can be slightly different. Consult with your doctor on determining what blood pressure is normal for you.

Before taking blood pressure, rest for approximately five minutes and take your blood pressure while relaxed in a quiet environment.

Measure blood pressure using the correct posture and do not move or speak during measurement.

Avoid exercise, eating, drinking alcohol, smoking and other activities that affect your blood pressure right before a measurement.

Take your blood pressure at the same time every day.

The ambient temperature should be approximately 20°C when you take your blood pressure.

Reading may vary slightly, depending on posture, whether sitting or lying.


## PREPARATION BEFORE USE

Install batteries and activate the clock before use. Measurement can be made with clock being inactivated. However, activate the clock to save measurement date and time with readings for better review. The clock is inactivated as the batteries are removed.

The product can be also operated with AC adaptor. See page 8.

### 1 INSTALLING BATTERIES

Open the battery compartment cover, pressing the projection between [ PUSH ] and [  ] at the bottom of the main unit.

Install four "AA" type batteries into the compartment. Make sure that the polarities correspond to the (+) and (-) marks inside the battery compartment.

Batteries can be easily installed or removed as their (-) ends are pushed against the spring.


You may use nickel hydride rechargeable batteries.

NOTE: These batteries cannot be charged with this product. Use the designated battery charger. Also, for details regarding use of rechargeable batteries, follow the directions accompanying them.


You will hear a beep as the batteries are installed.


NOTE: The beep will not sound if AC adaptor is connected.

AA (LR6)  
BATTERIES

Close the battery compartment cover.

Do not force the battery cover into position. Insert two tabs to the main unit first, Fig.1). Then close the cover, pressing [  ] section Fig.2). Close the cover completely Fig.3).

Battery symbol [  ] will flash when battery power is low. Measurement can not be made while the symbol is displayed. Replace the batteries. All four batteries need to be replaced with new ones; do not mix new and old batteries or different batteries.

The symbol might appear only during measurement although you can review readings in memory. This is because more battery power is required for conducting a measurement, e.g. inflating the cuff, than displaying the readings in memory. Please have the batteries replaced.

*The enclosed batteries are for monitoring, and their life may be shorter than that of commercial batteries.*


## 2 ADJUST DATE AND TIME AND ACTIVATE THE CLOCK

Touch "SET" KEY and keep your finger on the key until "20 11" starts flashing.

Clock is set in the order of year, month, day, hour and minute.


The flashing number increases with MEMORY KEY [1] and decreases with KEY [2]. The number will be fast-forwarded if you keep your finger on the key. Touching "SET" key will fix the number and the next article will flash. Touching "START/STOP" key will terminate the setting.

As the date and time are set, clock is activated. Display of the clock while the device is turned off indicates that the clock is activated.

### OPERATION WITH AC ADAPTOR

AC CONNECTOR is located at the side of the monitor.

Use the designated AC adaptor only.

It is suggested to take out the batteries when using the adapter. However, the clock is inactivated as the AC adaptor is unplugged and therefore you may insert the batteries to keep the clock activated. In such case, please make a few measurement without plugging the adaptor to avoid battery leak.

Do not connect the AC adaptor to any power socket where accessibility is not satisfying.


Protection: This device is double insulated and protected against short circuit and overload by a primary thermal fuse.

Enclosures and Protective Covers: Equipment enclosed to protect against contact with live parts, and with parts which can become live (finger, pin, hook test).

## MAKING A MEASUREMENT

Insert AIR PLUG into AIR CONNECTOR before starting a blood pressure measurement.

AIR CONNECTOR

The cuff is applicable to arm with circumference between 22 and 42 cm. Take your arm size before measurement.

### 1 SIT ON A CHAIR AND APPLY CUFF.

Find a chair and a table so that you can be comfortably seated with your feet flat on the floor

and your back and arm supported.


The cuff condition cannot be determined properly unless the cuff is applied at the correct site in the correct direction.

Inaccurate readings may result if cuff is wrapped over heavy clothing or tight shirt that may restrict circulation in your upper arm or if sleeve is rolled up.


Put the cuff on your left arm with AIR HOSE positioned toward your hand.

Attach the cuff with hook and loop fastener with the lower edge of the cuff approximately 2 to 3 cm above the inner elbow.


The cuff should fit snugly, leaving enough space for two fingers to be put between your arm and the cuff.

If the cuff is wrapped more tightly or loosely than this, inaccurate blood pressure readings may result.

Adjust the cuff position to make sure that AIR HOSE comes out to the center of the forearm towards the palm.

Press the surface of the cuff to make sure that the hook and loop fastens securely.

**2** Let your forearm rest naturally on the table and keep the cuff at the height of your heart during measurement.

If the cuff is lower (higher) than the heart, the measured reading tends to become higher (lower).

Take deep breaths and relax.

Do not move or speak while taking your blood pressure.

Do not cross your legs during measurement.

Do not hold your breath.


**3** TURN ON THE MONITOR.

Touch "START/STOP" KEY.

Air is exhausted from the cuff.

Automatic inflation starts.

CUFF SYMBOL is displayed.

Touch "START/STOP" KEY again to cancel measurement.

The monitor will exhaust air from the cuff and turn off.

Inflation stops at optimal pressure and pressure starts to decrease.


PULSE SYMBOL is displayed as pulse is detected.

Air is released from the cuff and measured values are displayed as measurement completes.

Refer to page 3 for [  ].

Refer to page 4 for [  ].

Refer to page 12 for "Err" display.

#### 4 SELECT MEMORY BANK TO SAVE THE MEASUREMENT RESULT.

MEMORY SYMBOL is displayed below MEMORY KEY.

Touching MEMORY KEY will switch the bank.

MEMORY MARK

Readings are automatically saved in the bank, either MEMORY [1] or [2], selected when the monitor is turned off after measurement. Make sure to confirm the bank selected after a measurement.

Refer to page 4 for memory function.

#### 5 TURN OFF THE MONITOR.

Touch "START/STOP" KEY.

The monitor will be turned off. Even if you do not turn off the monitor, it turns off automatically after 3 minutes.

#### AUTOMATIC REPEATED INFLATION

If you move or strain your arm or hand, inflation may be repeated multiple times to obtain a measurement. Repeated inflation does not indicate malfunction of the monitor.

**Do not execute repeated measurements for congestion of blood could result in false measurement. Let your arm rest for at least 5 minutes.**

# VIEWING THE SAVED READINGS (MEMORY FUNCTION)

Turn off the monitor once with "START/STOP" KEY when the results are displayed after a measurement.


Refer to page 4 for memory function.

## 1 TOUCH MEMORY KEY.

Touch MEMORY KEY [1] to review readings saved in the memory bank [1] and MEMORY KEY [2] to review those in the bank [2].

The average of the saved readings is displayed.


Note: Average will not be displayed unless there are two or more readings saved.

## 2 TOUCH THE SAME MEMORY KEY AGAIN.

The latest reading saved in the selected bank is displayed.

The display at the top shows memory number, measurement date and time alternately. The latest reading is displayed as memory "No. 1".

The results are recalled from the latest to the oldest with operation of MEMORY KEY.

The bigger memory number indicates the older number.

## SWITCHING MEMORY BANKS

Touching MEMORY KEY [2] at display of reading in the bank [1] will switch the display to reading in bank [2] and touching MEMORY KEY [1] again will return to the display of result in the bank [1].

## 3 TURN OFF THE MONITOR.

Touch "START/STOP" KEY.

Even if you do not turn off the monitor, it displays the result for approximately 30 seconds and will turn off automatically.

## DELETING SAVED READINGS

### Deleting individual readings

Touch the memory key to display the latest set of readings. Then scroll through the readings by repeatedly touching the memory key until you reach the reading to be deleted. Then hold down memory key for 4 seconds until display flashes and continue to hold down key for about 4 seconds until reading is deleted and the display is cleared [ - - - ].

### Deleting all the readings in a memory bank.

Select and display the average readings for the memory bank to be deleted. Then delete by holding down the memory key continuously until the display is cleared [ - - - ] as described in the sentence above for the deletion of individual readings.

# TROUBLESHOOTING

???

Blood pressure is extremely high or low.

Measurement is not taken with correct posture or the cuff is wrapped incorrectly. **Reconfirm the measurement procedures.**

The reading was effected by moving or talking. **Remain still and quiet during measurement.**

Measurement was taken right after exercise or other activities influencing the reading. Measure again after resting for more than 5 minutes. **See TIPS ON TAKING YOUR BLOOD PRESSURE.**

???

Measured values vary all the time.

Mental and physical condition or measurement condition is influencing blood pressure. **Make measurement under the same conditions.**

???

Measured value is different from that taken at a hospital.

Mental condition such as nervousness at a hospital or relaxed feeling at home is influencing blood pressure. **Take records of blood pressure taken at home and consult with your doctor.**

Blood pressure is changing 24 hours a day in relation with mental and physical conditions. Eating, drinking, smoking, exercising or bathing will alter blood pressure. Blood pressure will also be influenced by mental condition such as feeling tense or relieved and health condition. It is important to know the trend of blood pressure by making measurement regularly. Make your own rules including conditions and time to make measurement and take blood pressure everyday.

???

Inflation is repeated.

The cuff is inflated again if sudden change of pressure, which could be due to movement or straining of arm or wrist, was detected. Inflation is repeated until blood pressure is taken. **Repeated inflation is not a malfunction of the monitor. Do not move or speak during measurement.**

Err  
300

[Err 300] is displayed.

Over pressure: blood pressure could not be taken due to moving or speaking although the cuff was inflated to the maximum pressure. Do not move or speak during measurement.

Err  
-1

Err  
-3

[Err -1] or [Err -3] is displayed.

Blood pressure could not be taken due to moving or speaking. **Do not move or speak during measurement.**

Err  
-2

[Err -2] is displayed.

The cuff is not securely connected to the monitor. **Secure the air connector.**

The cuff is not wrapped appropriately. **Reconfirm how to put the cuff on your arm.**

[] is displayed.

Batteries are weak. **Replace all batteries with new ones.**

0  
0

[0 0] is displayed.

"START/STOP" KEY was being touched when batteries are being installed or when AC adaptor is being plugged in. **Turn off the monitor once with "START/STOP" KEY.**

???

Nothing is displayed.

Batteries are depleted. **Replace all batteries with new ones.**

Batteries are inserted incorrectly. **Reinsert the batteries correctly.**

Battery terminals are not clean. **Clean the terminals with dry cloth.**

AC adaptor is not connected. **Connect AC adaptor.**

Touching the keys with tip of finger or nail is touching the key. **Touch the key with the plane of your finger.**

???

Clock is not displayed.

The clock is inactivated. NOTE: Clock is inactivated when batteries and/or AC adaptor are/is removed. **Adjust date and time and activate the clock.**

???

Measurement date and time are displayed with [ -/-/ - ] and [ -:-: - ].

The clock is not activated. **Adjust date and time and activate the clock.**

Measurement was taken before the clock was activated. **Measurement date and time cannot be saved without clock being activated.**

If you cannot get correct measurement with the methods above, contact your dealer. Do not disassemble or tamper with the internal mechanism.

## TECHNICAL DESCRIPTION FOR ELECTROMAGNETIC DISTURBANCES

DSK-1031 complies with the Electromagnetic Disturbances standard, IEC60601-1-2:2014. As a medical electrical equipment, special precautions regarding the electromagnetic disturbances shall be taken at usage of the device according to the information provided below.

- The device is not intended for use in environments where the intensity of electromagnetic disturbance is high, such as near active HF surgical equipment and MRI (magnetic resonance imaging) equipment etc.
- Use of the device adjacent to or stacked with other equipment must be avoided because it could result in improper operation.

- Use of accessories other than those specified or provided by the manufacturer could result in increased electromagnetic emissions or decreased electromagnetic immunity of the device and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used at least 30cm away from any part of the device, including specified cables. Otherwise, degradation of the performance of this equipment could result.

Please contact your dealer or the manufacturer for specific information regarding the compliance to the standard.

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#### WARRANTY

This equipment is guaranteed for the period of 2 years after the date of purchase against manufacturing defects when returned along with the proof of date of purchase to the dealer from whom the purchase was made. During this period, the unit will be repaired or replaced free of charge if the failure is attributable to faulty design or manufacture. This warranty does not cover damage or malfunctions caused by improper handling or use contrary to the instructions in this manual. Please contact your dealer for additional warranty provisions which may remain effective after the manufacturer's warranty period has expired.

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**NIHON SEIMITSU SOKKI CO., LTD.**
**MDSS GmbH**