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
Instruction Manual NISSEI

Operating Principles

The cuff applies pressure to the artery in your arm which pulses with your heart rate. The pressure in the cuff moves with the pulse. This pulse changes according to the amount of pressure the cuff applies to your arm. The oscillometric method can find the systolic and diastolic values as the level of pressure applied from the cuff changes gradually.

© 2014 NIHON SEIMITSU SOKKI CO., LTD.
Reproduction and reprinting are prohibited.

PRECAUTIONS AND CONTRAINDICATIONS

Read this instruction manual carefully before use and use this product correctly. The instructions below are to help you use the product correctly and to prevent you from being harmed or injured.

- Prohibition** The subjects or use which exceeds the responsibility of the manufacturer such as the use beyond the design limit or the unjustified use. (Ones with possible death or serious injury are indicated.)
- Caution** General cautions for use of the product. (The subjects which may result in personal injury or property damage in case the product is incorrectly handled are indicated.)

*Property damage means consequential damage for house, household belongings, livestock or pet.

Prohibition


- Do not use the product for infants or persons who cannot express their will. It may cause an accident or an injury.
- Do not bring the product in the test room when having a MRI scan. It may cause the attraction by magnetic force or the thermal burn.
- Do not use the product near flammable gas, such as anesthetic gas, or in locations or atmospheres with high oxygen concentration such as hyperbaric oxygen chambers/therapy devices or oxygen tents. It may cause malfunction, breakage, deterioration or induction of explosion.
- Do not make any diagnoses or treat yourself based on your own measurement results. Make measurement based on the consultation with your doctor and receive professional treatment and prescription.
- Do not use the product, if you have an injury on/around your arm. It may cause the aggravation of the injury.
- Do not use the product during intravenous drip injections or blood transfusions. It may cause an accident or an injury.
- Do not use the product with exceeding its service life. The service life of the product is either 5 years or 30,000 times.
- Do not use the product in the medical institution or public area where the product is subject to be used by unspecified persons. It may cause an accident or a trouble.

Caution

- Do not use the product near mobile phones, micro ovens or other strong electromagnetic devices. It may cause malfunctions or damage the product.
- Be cautious about using the product in conjunction with the medical electrical equipment listed below. It may cause malfunctions.
 - Pace makers, implantable defibrillators, or other wearable medical electrical equipment that is easily affected by electromagnetic force
 - Wearable medical electrical equipment such as an electrocardiographic monitor.
- Consult a doctor when persons in pregnancy use the product.
- Do not use the product for purposes other than measuring blood pressure. It may damage the product or cause an accident.
- Do not use the cuff around anywhere other than your upper arm. It may damage the product or cause an accident.
- When using the monitor in an environment where systems connecting fluid with blood vessels such as dialysis and intravenous drip injection are used, do not connect the blood pressure monitor and the system. An accident may be caused if the air enters a blood vessel.
- Consult a doctor before using the product for persons suffering from arrhythmia or peripheral circulatory disorders caused by diabetes, liver disease, arteriosclerosis or hypertension.
- When using the cuff, use the supplied equipment. Using equipment other than that supplied may cause measurement error.
- Measure the circumference of your upper arm and make sure that it is within the suggested range. Using the product for an arm outside the suggested range may cause measurement error.
- Do not place heavy items on the product or tie up the air hose. It may damage the product.
- Do not make an impact on or drop the product. It may cause malfunctions or damage the product.
- Wear the cuff correctly and keep it level with the heart during measurement. Wearing the cuff incorrectly or at the wrong height may cause measurement error.
- Be cautious that blood pressure is affected by the factors below.
 - Time and season
 - Physical activity
 - Taking medicine for high blood pressure
 - Eating and drinking (including Alcohol)
 - Smoking
 - Mental nervousness
 - Bathing
 - The urge to urinate
 - Conversation
 - Other environment (during medical examination etc.)
 - Measurement posture (including the cuff level against the heart)
- Press the START/STOP BUTTON if the product applies an excessive level of pressure to your arm. Pressing the START/STOP BUTTON releases the air from the cuff immediately, and the product stops measuring.
- The cuff may leave bruises on your arm. If you feel any pain, press the START/STOP BUTTON to stop the operation of the product. Pressing the START/STOP BUTTON releases the air from the cuff immediately, and the product stops measuring.
- Do not plug/unplug the AC adaptor with wet hands, or touch it with tweezers, or other electronically conductive devices.
- Do not touch another person while touching the batteries, battery terminals of the device or the output plug of the AC adaptor, such as when replacing the batteries or plugging/unplugging the AC adaptor.
- Do not use new batteries in conjunction with old batteries or different types of batteries at the same time. It may generate heat and damage the product.
- Pull the plug out from the power outlet if anything unusual is detected during use. Plugging the AC adaptor in a power outlet located behind furniture may be inconvenient if it needs to be pulled out in an emergency.
- Do not unplug the AC adaptor by pulling the cable. It may damage the cable and cause a fire or electric shock.
- Take the batteries out when not using the product for a long period of time. Battery fluid may leak from the batteries and damage the product.
- Do not allow children to handle the product alone or leave the product in a place children can reach. It may injure the children.
- Do not disassemble, repair or modify the product. It may cause malfunctions or damage the product.
- Check if the product is dirty or wet before use.
- If the product stops functioning properly, stop using it. Please contact us for repair or assistance.
- Follow the local regulations of your area when disposing the product or batteries.

PART NAMES AND PRODUCT COMPONENTS

Make sure that the following items are in the package. If you do not have all the necessary items in the package, please contact the shop from which you brought the product or call our customer service (☎ 0120-211-164).

INSERTING BATTERIES

- Remove the cover of the battery compartment on the underside of the unit.
- Insert four AA type batteries into the compartment, observing proper polarity (refer to the marking in battery compartment). Inserting the battery with (+) first may damage the spring.
- Close the battery compartment by pressing the battery cover back into the unit until it clicks.

The battery symbol (🔋) flashes when the batteries are low. The main unit cannot be used if the battery symbol is lit (not flashing). Replace all the batteries with new AA batteries when the battery symbol is displayed.

The product can also be used with the optional AC adaptor. Connect the AC adaptor to the AC ADAPTOR JACK located to the right of the monitor. The AC Adaptor is sold separately. Purchase the designated AC adaptor, ADP-W5, separately from the unit.

ADJUST DATE AND TIME AND ACTIVATE THE CLOCK

When the batteries are properly inserted, "20 14" starts flashing on the display. Set the year, month, day, hour and minute in order.

- Set the year**
The number increases when the MEMORY BUTTON is pressed. Touch the START/STOP BUTTON to confirm the setting. Proceed to the next setting.
- Set the month, day, hour, and minute in the same way**
After adjusting the number, press the START/STOP BUTTON to confirm. When the minute is set, the procedure is complete.

If you wish to change the time, take the batteries out and press the START/STOP BUTTON. Make sure there is nothing shown on the display and then put the batteries back in.

TIPS ON TAKING YOUR BLOOD PRESSURE

The reading may vary slightly, depending on your posture. Take your blood pressure at the same time every day.

- Designate a time in the morning and the evening.**
It is said that your blood pressure is most stable before going to bed and soon after waking up in the morning. Check your blood pressure in the morning within an hour of waking up, before eating, and in the evening an hour or more after eating.
- Use one designated arm (right/left) when taking blood pressure.**
The blood pressure of one arm may differ from that of the other. Designate an arm for measurement when you take your blood pressure.
- Take your blood pressure in a quiet environment and use the same environment every time.**
Changes in the surroundings such as movement or noise may affect your blood pressure. Your blood pressure may increase when the temperature is low. A room temperature of 20° is ideal.
- It is best to take your blood pressure when you are relaxed.**
Mental tension or nervousness may affect your blood pressure. Take a deep breath a few times, rest for approximately five minutes, and take your blood pressure.
- Sit quietly and stay still while taking your blood pressure.**
Your blood pressure can be affected by speaking and moving your body. This product uses the oscillometric method for measurement of blood pressure and pulse rate. The oscillometric method can detect subtle pulsatory changes. For this reason, for accurate results, it is best not to speak, or move your body while taking your blood pressure.

Note
If you need to urinate, do so before taking your blood pressure.
Do not take your blood pressure after smoking.
Rest after taking a bath or exercising before taking your blood pressure.
Your blood pressure increases if you need to urinate, and after eating, smoking, or exercising.

ATTACHING THE CUFF

- Connect the cuff and the main unit**
Insert the AIR PLUG into the AIR CONNECTOR securely.
- Sit on a chair and attach the cuff**
Put the cuff on your left arm with the AIR HOSE positioned toward your hand.
* If you are unable to wear the cuff on your left arm, wear it on your right arm. Once you have decided to wear it on one arm, always wear it on the same arm.
- Fasten the cuff**
Fasten the cuff with the velcro, making sure that the lower edge of the cuff is approximately 2 to 3 cm above the inner elbow. Hold the cuff around your arm and fasten it tightly.

If the cuff still moves around, fasten it more tightly. Be careful not to fasten it too tightly.

- Adjust the position of the cuff**
Adjust the position of the cuff so that the ARTERY MARK is over the brachial artery or the AIR HOSE is at the center of your arm.

- Let your forearm rest naturally on the table and keep the cuff at the height of your heart during measurement**
If the cuff is higher than the heart, the measured reading tends to be lower whereas if the cuff is lower than the heart, the measured reading tends to be higher.

In order to measure your blood pressure accurately, adjust the height of your chair or place a cushion under your elbow so that the cuff is level with your heart.

You may lay on your back and take your blood pressure. Even then, keep the cuff at the same level as your heart. The measured value may vary depending on your posture (whether sitting or lying on your back).

MAKING A MEASUREMENT

- Start measuring your blood pressure**
Press the START/STOP BUTTON. All the indicators light up and measuring begins automatically.
- The measured value displayed**
The cuff releases the air, the measurement stops and the main unit displays the result.
- Turn off the main unit**
Pressing the START/STOP BUTTON saves the measured value, and then the main unit turns off automatically.

- BODY MOTION SYMBOL**
When you move your body during measurement, the result may not be as accurate as it may be otherwise. Sit still and take your blood pressure again.
- IRREGULAR PULSE RHYTHM SYMBOL**
When the intervals between pulses are irregular, the IRREGULAR PULSE RHYTHM SYMBOL is displayed.

VIEWING AND DELETING THE SAVED READINGS

The main unit can save up to 60 readings and their average value.

- Display the readings**
Press the MEMORY BUTTON. The average blood pressure value is displayed. Press the MEMORY BUTTON repeatedly to see the other results from the latest to the oldest.

The result number, date and time of measurement change automatically.

- Delete the saved results**
1. Press the MEMORY BUTTON and select the results you wish to delete.
* Delete all the results by displaying the average.
2. Press and hold the MEMORY BUTTON down for more than 4 seconds.
3. Hold the MEMORY BUTTON down until the results start flashing.
4. The results are deleted.

You cannot restore deleted results.

TROUBLESHOOTING

When	Cause	Solution
Blood pressure is extremely high or low.	The cuff is not level with your heart.	Keep the cuff at the level of your heart.
	The cuff is wrapped incorrectly.	Make sure that the cuff is wrapped correctly.
	You moved or spoke during measurement.	Do not move or speak during measurement.
The results are different from usual.	The condition of your mind and body and the measurement method are affecting the result.	Take your blood pressure in the same environment as usual.
The measured value differs from one taken at a hospital.	Mental conditions such as nervousness at a hospital or feeling relaxed at home is influencing blood pressure.	Take records of your blood pressure taken at home and consult with your doctor.
Err is displayed.	You moved or spoke during measurement.	Do not move or speak during measurement.
Err2 is displayed.	The AIR HOSE is bent or squeezed.	Make sure that the AIR HOSE is not squeezed or bent.
	The AIR CONNECTOR is not connected properly.	Make sure that the AIR CONNECTOR is connected properly.
	The cuff is wrapped incorrectly.	Reconfirm how to put the cuff on your arm.
Err300 is displayed.	You moved or spoke during measurement.	Do not move or speak during measurement.
	The pressure reached its maximum limit as you moved your hand or arm, or spoke during measurement.	Do not move or speak during measurement.
Nothing is displayed even when the buttons are pressed. (Measurement does not start)	The batteries are weak.	Replace all the batteries with four new batteries.
	Batteries are inserted incorrectly.	Reinsert batteries correctly.
	The batteries terminals are not clean.	Clean the terminals with a dry cloth.
Err is displayed.	The buttons were being pressed when the batteries were put in or when the AC adaptor was plugged in.	Press the START/STOP BUTTON then turn off the main unit.

If the above solutions do not fix the problem, contact the store from which you bought the product.

GENERAL INFORMATION

Blood pressure
Blood pressure is the pressure exerted by circulation of blood upon the walls of blood vessels.


Systolic blood pressure value
When the heart shrinks and a larger amount of blood flows out, blood pressure increases.

Diastolic blood pressure value
When the heart expands and blood returns to the heart, blood pressure decreases.

Blood pressure is constantly changing
Blood pressure varies depending on both mental and physical conditions.
Blood pressure is affected by the factors below.
• Time and season
• Physical activity
• Taking medicine for high blood pressure
• Eating and drinking (including Alcohol)
• Smoking
• Mental nervousness
• Bathing
• Conversation
• The urge to urinate

In order to take your blood pressure accurately, rest longer than 5 minutes and measure it a few times.

High blood pressure
Measured blood pressure values are classified against WHO guideline. Measured blood pressure values vary depending on racial background, area, sex, and age. Therefore the ideal blood pressure is different for each individual. Do not diagnose yourself based on your own measurement results. Consult with your doctor.

	Systolic blood pressure	Diastolic blood pressure
Optimal	< 120	and < 80
Normal	120 - 129	and/or 80 - 84
High Normal	130 - 139	and/or 85 - 89
Grade1 Hypertension	140 - 159	and/or 90 - 99
Grade2 Hypertension	160 - 179	and/or 100 - 109
Grade3 Hypertension	≥180	and/or ≥110
Isolated Systolic Hypertension	≥140	and < 90

WHO/JSH Hypertension guidelines (2009)

CARE AND MAINTENANCE

This product can be used without the professional skills or knowledge required to perform maintenance on the product. The overall condition of the product may deteriorate depending on the number of times it is used, the way it is used/maintained and the number of years it is used. The service life of the product is either 5 years or 30,000 times.

- Cleaning the product**
Make sure that you clean the product after use. If it becomes dirty, wipe the surface of the main unit with a cloth using warm water or soapy water, then wipe it dry with a soft cloth. To clean the cuff, lightly pat the surface of it with a cloth using neutral detergent. Make sure that water never enters the AIR HOSE. Dry the cuff thoroughly. Do not use a washing machine.
Do not use benzene, gasoline, paint thinner or alcohol. It may damage the product.
- Storing the product**
Do not store the product in the places described below for a long period of time. It may damage the product or may allow it to deteriorate.
 - Places exposed to direct sunlight
 - Places affected by differences in temperature or warm places with a high degree of humidity
 - Dusty places
 - Wardrobes where insecticide is used

When not using the product for a long period of time, remove the batteries from the main unit. Otherwise, battery fluid may leak from the batteries and damage the product.

PRODUCT SPECIFICATIONS

Model : BLOOD PRESSURE MONITOR DS-B10
 Operating Principle : Oscillometric method
 Display : 13 digits liquid crystal display
 Pressure Indication Range : 3 to 300 mmHg (cuff pressure)
 Error of cuff pressure indication : within ±3mmHg
 Measuring Range : 50 to 250 mmHg (systolic)
 40 to 180 mmHg (diastolic)
 40 to 160 bpm (pulse rate)

Errors in blood pressure in clinical testing*
 • Auscultation : Average Error within ±5 mmHg
 Standard deviation : less than 8 mm

Pulse rate measurement accuracy : within 5%
 Pressure : Automatic inflation
 Exhaust : Automatic exhaust valve
 Power Supply : 4 AA batteries (alkaline) or AC adaptor (ADP-W5)
 Main Unit Weight : Approx. 220 g (without batteries)
 Power Consumption : When using 4 AA batteries (alkaline): 6 V DC 4 W
 When using AC adaptor ADP-W5: Power consumption for AC adaptor ADP-W5 100-240 V AC 50/60 Hz 0.12 A
 Power consumption for main unit 6 V DC 500 mA

Memory : The main unit stores 60 readings and their average
 Operating Conditions : +10°C to +40°C, Relative humidity 30% to 85% RH (no condensation)
 Transportation/Storage Conditions : -20°C to +60°C, Relative humidity 10% to 95% RH (no condensation)
 Cuff coverage arm circumference : Approx. 22 to 42 cm
 Main Unit Dimensions : Approx. 122.2 x 116.5 x 69.9 mm (W x D x H)
 Main Unit Weight : Approx. 220 g (without batteries)
 Protection against electric shock : Internally powered equipment/Class II equipment
 Protection against harmful ingress of water : IPX0
 Classification of safety level when used in air/flammable anesthetic gas, oxygen, nitrous oxide/flammable anesthetic gas.
 • Not suitable for use in these environments.
 Classification by function (operation) mode
 • Continuous function (operation) equipment
 • Class II equipment
 • GF Applied Parts
 • Read the attached manual.

Specifications are subject to change without notice due to improvements in performance.
 This product conforms with IEC60601-1-2:2007 (Classification and CISR classification group: Group 1, Class B). Please contact the manufacturer for detailed information.
 * This measuring method is based on ANSI/AAMI SP10:2002, American National Standard for electronic or automated Sphygmomanometers, and Amendment.