NISSEI

**DIGITAL
BLOOD PRESSURE
MONITOR** *Model DS-1902*

**INSTRUCTIONS
ENGLISH**

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CONTENTS

PARTS NAMES AND PRODUCT COMPONENTS	4
PICTORIALS	5
BASIC PRODUCT INFORMATION	11
PRECAUTIONS FOR USE	12
TIPS ON TAKING YOUR BLOOD PRESSURE	12
INSTALLING BATTERIES AND SETTING THE CLOCK	13

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
PARTS NAMES AND PRODUCT COMPONENTS

- | | | | | |
|------------------|------------------|----------------------|----------------------|----------------------|
| ① MONITOR | ② CUFF | ③ 4 AA BATTERIES | ④ CARRYING BAG | ⑤ INSTRUCTION MANUAL |
| ① MONITOR | ② CUFF | ③ 4 AA BATTERIES | ④ CARRYING BAG | ⑤ INSTRUCTION MANUAL |
| ② CUFF | ③ 4 AA BATTERIES | ④ CARRYING BAG | ⑤ INSTRUCTION MANUAL | |
| ③ 4 AA BATTERIES | ④ CARRYING BAG | ⑤ INSTRUCTION MANUAL | | |
| | | | | |

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
PULSE RATE
PULSRATE **POULS**
BATTITO CARDIACO
RAPIDEZ DEL PULSO
WARTOŚĆ PULSU

/min


IRREGULAR PULSE RHYTHM
INDICATION

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








BASIC PRODUCT INFORMATION

This manual is intended to assist the user in the safe and efficient operation of Digital Blood Pressure Monitor DS-1902. The product must be used in accordance with the procedures contained in this manual and must not be used for purposes other than those described herein. It is important to read and understand the entire manual. In particular, please read carefully and become familiar with the section entitled "TIPS ON TAKING YOUR BLOOD PRESSURE".

Indications for Use

The DS-1902 system is intended for the non-invasive measurement of systolic and diastolic blood pressure and determination of pulse rate in adults, i.e., age 12 and above; this unit is not designed for neonatal use. Also, an inaccurate reading may result if this instrument is used on a child's arm. Consult your physician if you wish to take a child's blood pressure. The product is recommended for use by patients with labile blood pressure or known hypertension in a home care environment as an adjunct to medical management. The cuff will accommodate an upper arm circumference range of approximately 23 to 43 centimeters. Pressure is measured over a range of 0 to 300 mmHg and pulse rate over a range of 40 to 160 beats/minute.

Method of Measurement

This product employs the oscillometric method for measurement of blood pressure and pulse rate. The cuff is connected to the main unit and wrapped around the arm. As START/STOP BUTTON is pressed, the system starts automatic pressurization, during which measurement of blood pressure is done. Circuits within the cuff sense the small oscillations in pressure against the cuff produced by the expansion and contraction of the arteries in the arm in response to each heart beat. Pressurization stops when the cuff is inflated enough to determine diastolic and systolic blood pressures and air is exhausted from the cuff. The amplitude of each pressure waves is measured, converted to millimeters of mercury, and displayed on the LCD as a digital value. If the pulse rhythm detected during recent measurement was irregular, irregular pulse rate rhythm indication be displayed. A memory circuit stores the 30 most recent measurement results with date and time for comparison and computes the average value of stored data.

DS-1902 meets the stringent requirements of German Hypertension League to carry its quality marking.

PRODUCT SPECIFICATIONS

Model	: DS-1902	Operating Environment	: +10°C to +40°C, 85% relative humidity or below
Operating Principle	: Oscillometric method	Storage Environment	: -5°C to +50°C, 85% relative humidity or below
Indicator	: 12 digits liquid crystal display	Cuff	: Coverage arm circumference ; 23 to 43cm
Pressure Indicating Range	: 0 to 300 mmHg (Cuff pressure)	Main Unit	: Weight ; Approx. 330 g (without batteries) Size ; 148 x 115 x 54 mm (W x D x H)
Measuring Range	: 50 to 250 mmHg (systolic) 40 to 140 mmHg (diastolic) 40 to 160 bpm (pulse rate)	Key to symbols	
Accuracy	: ± 3 mmHg (Cuff pressure), $\pm 5\%$ of reading (pulse rate)		: Type BF equipment
Inflation	: Automatic inflation		: Important ; Read operating instructions
Exhaust	: Automatic exhaust valve		
Power Supply	: 4 pcs. 1.5 volt "AA" (R6P/LR6) type batteries or AC adapter		
Power Consumption	: 4W (max.)		
Memory	: 30 measurements x2 and averages		
			This device complies with EMC (IEC60601-1-2), EN1060-1 and EN1060-3. Specifications are subject to change without notice due to improvements in performance.

PRECAUTIONS FOR USE

Do not use cuffs other than the original cuff included with this product.

The system may fail to yield specified measurement accuracy if operated or stored in temperature or humidity conditions outside the limits stated in the specifications section of this manual.

Do not use this instrument without consultation with your doctor if you are under dialysis therapy or on anticoagulants, antiplatelets or steroids. Use of this instrument under such conditions could cause internal bleeding.

To avoid any possibility of accidental strangulation, keep this unit away from children and do not drape AIR HOSE around your neck.

For specific information on your blood pressure, contact your physician. Never make any judgment on your own regarding measurement results.

Use of this device in areas near portable phones, microwave ovens or other devices with strong electromagnetic field may cause malfunctions.

TIPS ON TAKING YOUR BLOOD PRESSURE

Blood pressure is a measurement of the force exerted by the heart in pumping the blood through the arteries and the resistance by the veins in this flow.

Blood pressure varies all the time, influenced by mental and physical factors and is never constant.


In general, blood pressure is highest during the working hours and gradually decreases during the afternoon and evening hours. It is low during sleep and increases at a relatively fast rate after arising from bed.

Causes for Changes in Blood Pressure

- Body movement
- Emotions
- Recent Urination or Bowel Movement
- Changes in the surroundings such as movement or noise, etc.
- Conversation
- Eating
- Drinking Alcohol

- Mental Tension
- Smoking
- Temperature

Blood pressure changes shown below.
Upper arm blood pressure data

Before taking blood pressure, rest for approximately five minutes and take your blood pressure while relaxed in a quiet environment.

Measure blood pressure using the correct posture and do not move or speak during measurement.

Avoid exercise, eating, drinking alcohol, smoking and other activities that affect your blood pressure right before a measurement.

Take your blood pressure at the same time every day.

The ambient temperature should be approximately 20°C when you take your blood pressure.

Blood pressure measured at home tends to be lower than when measured in a hospital, clinic or doctor's office.

This is because you are tense at the hospital and relaxed at home. It is important to know your stable normal blood pressure at home.

Let a qualified physician interpret your blood pressure readings.

Depending on your age, weight and general condition, blood pressure can be slightly different. Consult with your doctor on determining what blood pressure is normal for you.

180~	(severe)					Hypertension
160~179	(moderate)					
140~159	(mild)					
130~139	High Normal					
120~129	Normal					
~120	Optimal					
SYS mmHg	80	85	90	100	110	
DIA	80	84	89	99	109	

Blood pressure classification – WHO (1999)

INSTALLING BATTERIES AND SETTING THE CLOCK

fig.1

1. Open the battery compartment cover.
2. Install four “AA” type batteries into the compartment.
Make sure that the polarities correspond to the (+) and (-) marks inside the battery compartment.

The batteries can be easily installed by pressing their (-) side against the spring.

3. Close the battery compartment cover.

Do not force the battery cover into position.

BATTERY REPLACEMENT INDICATION

Replace all the batteries when the battery replacement indication appears on the display or nothing is displayed.

Use alkaline batteries to increase battery life. Use of regular Carbon-Zinc batteries will necessitate more frequent replacement. *The enclosed batteries are for monitoring, and their life may be shorter than that of commercial batteries.*

4. Year flashes on the display.

Clock can be set only after the batteries are reinstalled. Adjust the clock so that measurement results are stored with correct date and time.

Use **MEMORY BUTTON** to increase the flashing year and **MEMORY SELECTOR** to decrease the year.

Press **START/STOP BUTTON** to confirm and move to next step.

5. Month flashes.

Adjust the month with **MEMORY BUTTON** and/or **MEMORY SELECTOR**, press **START/STOP BUTTON** to confirm.

6. Day flashes.

Adjust the day with **MEMORY BUTTON** and/or **MEMORY SELECTOR**, press **START/STOP BUTTON** to confirm.

7. Hour flashes.

The clock is in 24 hour mode.

Adjust the hour with **MEMORY BUTTON** and/or **MEMORY SELECTOR**, press **START/STOP BUTTON** to confirm.

8. Minute flashes.

Adjust the minute with **MEMORY BUTTON** and/or **MEMORY SELECTOR**, press **START/STOP BUTTON** to confirm.

Clock is set and the unit is turned off.

The clock is displayed while the unit is turned off.

The used electrical and electronic products are not household waste. Follow your national/local recycling rules to dispose of them properly. In the EU countries, please refer to waste management symbol(s) marked on the package or the instrument.

Operation with an AC adapter

AC ADAPTER JACK is located on the right side of the unit.

Do not use any adapters other than the one specified below.

Model	: JPI2245
Output voltage	: 6 VDC ±5%
Input voltage range	: Standard voltage ±15%
Maximum electric load	: 600mA or higher
Output plug polarity	: Center Minus

Remove the batteries when operating the unit with AC adapter for extended period of time.

Clock is reset when AC adapter is unplugged.

CORRECT MEASURING POSTURE

fig.2

Sit at the table and let the table support your arm as you take the blood pressure measurement.

Make sure that the measurement location on the upper arm is at approximately the same height as the heart, and that the forearm is extended naturally on the table and does not move.

You may lie on your back and take the measurement. Look at the ceiling, stay calm, and do not move your neck or body during the measurement. Again, make sure that the measurement location on the arm is at approximately the same height as the heart.

Measured data may vary slightly depending on the posture during measurement.

If the cuff is lower (higher) than the heart, the measured reading tends to become larger (smaller).

WRAPPING THE CUFF

fig.3

- 1. Place the cuff on your left arm with the air hose positioned toward your hand.**
- 2. Wrap the cuff around your arm with the edge of the cuff approximately 2~3 cm above the elbow. AIR HOSE should be on the palm side of the arm.**
- 3. Press the surface of the cuff to make sure that the hook & loop fastens securely.**

When wrapping the cuff, wrap it loosely enough around the arm so that two fingers can be placed between the cuff and the arm. If the cuff is wrapped more tightly or loosely than this, inaccurate blood pressure readings may result.

If you are wearing a shirt that might restrict circulation in your upper arm or you roll your sleeve up over the upper arm, the blood flow will be restricted, preventing accurate measurement.

MEASUREMENT PROCEDURES

fig.4

- 1. Insert TUBE PLUG into AIR CONNECTOR.**
- 2. Press START/STOP BUTTON.**

Deflation mark flashes.

Beep sounds when the cuff is deflated completely.

Automatic inflation starts and inflation mark flashes.

Inflation mark disappears and measurement begins.

Heart mark flashes with synchronizing beep as pulse is detected.

Long beep sounds as the measurement is complete.

Blood pressures and pulse rate are displayed.

The unit automatically exhausts the air from the cuff.

Take deep breaths and relax.

Do not move, chat or strain your arm or hand during measurement.

To stop measurement

Press START/STOP BUTTON and inflation stops, air is exhausted, and monitor turns off.

3. Press MEMORY SELECTOR and select M1 or M2 to save the result or not to save the result.

The result is stored in the bank selected when the unit is turned off. Irregular pulse rhythm indication will not be stored, only blood pressure and pulse rate are stored. When a measurement resulted in an error, it will not be stored.

4. Press START/STOP BUTTON to turn off the power.

If you forget to turn off the unit, it will be automatically turned off after 3 minutes.

Automatic Repressurization

If noise was detected or diastolic blood pressure suspected to be low, measurement is interrupted and cuff is deflated once and inflation starts again.

Do not execute repeated measurements for congestion of blood could result in false measurement. Let your arm rest for at least 5 minutes.

IRREGULAR PULSE RHYTHM INDICATION


A flashing heart mark in the measurement result display indicates irregular pulse rhythm. Irregular pulse rhythm can be a result of body movement during measurement or arrhythmias. Although continuous appearance of the indication under quiet measurements may suggest arrhythmias, do not make any judgment on your own before consulting with your doctor.

An example of irregular pulse rhythm compared to regular pulse rhythm


Measurement during inflation

Initially the cuff is rapidly pressurized to approximately 25 mmHg at the beginning of a measurement, during which the inflation mark is displayed. After the pressure reaches this point, the inflation mark disappears and the cuff is pressurized more gradually. During this gradual inflation, pulse is detected and heart mark flashes. Measurement ends when systolic blood pressure is determined and air is exhausted from the cuff.

"Personal Mode"

This blood pressure monitor employs what is called "Personal Mode" to shorten the measurement time and improves comfort. With "Personal Mode", the pressure to which the monitor rapidly inflates is adjusted based on the last three measurements. For example, high diastolic blood pressure values from past three measurements will lead to higher inflated pressure value than the initial value. In the event that this value is too high for your current blood pressure, the cuff will rapidly deflate and measurement will restart.


"Personal Mode"

To use "Personal Mode" efficiently, select your memory bank before starting a measurement; press MEMORY SELECTOR to see which memory bank is selected, press MEMORY SELECTOR button to change the memory bank number, then press START/STOP BUTTON once to turn off the monitor, and then press START/STOP BUTTON to start a measurement. Select  to avoid using the "Personal Mode".

MEMORY FUNCTION*fig.5*

After a measurement, to add current reading to memory, press START/STOP BUTTON once while measurement result is being displayed. Measured result is stored in either M1 or M2. Make sure desired memory bank is selected before turning the monitor off after a measurement. The result will not be saved when  is selected. Each of two memory banks, M1 and M2, can hold up to 30 results and their average. When the number of measurements exceeds 30, the oldest data will be deleted to record new data.

**1. Press MEMORY BUTTON to see stored data.**

Press START/STOP BUTTON once while measurement result is being displayed after a measurement.

The average of the stored results is displayed.

When there is only a single result, this reading will be displayed as the average.

The memory data is displayed for approximately 30 seconds. Approximately 30 seconds after MEMORY BUTTON is released, the apparatus will turn off.

2. Press MEMORY SELECTOR to change memory bank.**3. Press MEMORY BUTTON to move to next data.**

The memory data number 1 is the latest among the stored data.

Indication at the top alternates from memory data number to date and to time of measurement.

Every depression of the key switches among the memory data.





As MEMORY BUTTON is pressed, the memory data number increases; the bigger the number, the older the result.

4. Press START/STOP BUTTON to turn off the unit.**DELETING STORED DATA****1. Display the data to be erased or the average to erase all the data in the memory.****2. Press and hold down MEMORY BUTTON.**

The displayed data starts flashing.

3. Hold down the button until measurement result disappears.

ERROR DISPLAYS AND TROUBLESHOOTING

ERROR SYMBOL/SYMPTOM	CAUSE / CHECK POINT	REMEDY
	<p>OVER-PRESSURIZATION; The cuff was inflated to the maximum pressure because of movement of body etc.</p>	<p>Do not move during measurement.</p>
	<p>MEASUREMENT ERROR; Measurement could not be made because of moving or talking during measurement.</p>	<p>Remain still and quiet during measurement.</p>
	<p>INFLATION ERROR; TUBE PLUG is not correctly inserted. Noise was detected.</p>	<p>Reinsert TUBE PLUG and make sure that it is securely inserted. Do not move during measurement.</p>
	<p>NO POWER; Batteries are exhausted.</p> <p>Have the batteries polarities been positioned incorrectly? Are the battery terminals clean? AC adapter is not properly connected.</p>	<p>Replace all batteries with new ones. Reinsert the batteries in the correct position. Clean with a dry cloth. Check the AC adapter connection.</p>
	<p>LOW BATTERY; Battery is weak.</p>	<p>Replace all batteries with new ones.</p>
	<p>IMPROPER OPERATION; START/STOP BUTTON was accidentally pressed during battery replacement.</p>	<p>Press START/STOP BUTTON to turn off the power once and press it again to start a measurement.</p>

ERROR SYMBOL/SYMPTOM	CAUSE / CHECK POINT	REMEDY
Measurement is interrupted once and cuff is deflated and inflated again.	When noise is detected or diastolic blood pressure is low, the cuff is deflated and inflated again. Did you move during measurement?	This does not indicate a problem. Do not move during measurement.
Cannot make great number of measurements.	Are you using alkaline batteries?	Use alkaline batteries to increase battery life.
Cannot complete the measurement.	Have the batteries run down?	Replace all four batteries with new ones.
Blood pressure is different each time. The reading is extremely low (or high).	Are you measuring with correct posture? Blood pressure readings constantly vary with time of measurement and nervous condition.	Measure with the correct posture. See TIPS ON TAKING YOUR BLOOD PRESSURE.
Pulse rate is too low (or too high).	Did you move during the measurement? Did you take measurement immediately after exercise?	Remain still during the measurement. Measure again after resting for more than 5 minutes.
Measurement result is not stored.	Is it the right memory bank? Did you select correct memory bank number after the measurement?	Confirm the memory bank number. Make sure desired memory bank is selected before turning off the monitor after a measurement. When is  displayed the result will not be saved.
The power is automatically turned off.	Have you left the instrument untouched after the measurement?	This does not indicate a problem. The instrument automatically turns off 3 minutes after a measurement or 30 seconds after memory display.

If you cannot get correct measurement with the methods above, do not tamper with the internal mechanism. Contact your dealer. If the unit malfunctions, please return it to the dealer or an authorized service representative for service according to the warranty.

CARE AND MAINTENANCE

Because the unit includes precision parts, care should be taken to avoid extreme temperature variations, humidity, shock, dust, and direct sunlight. Do not drop or strike the unit. Make sure not to expose the unit to moisture. **This unit is not water resistant.**

Use only a soft, dry cloth to clean the unit. Do not use gasoline, paint thinner, or other strong solvents on the unit. Since the cuff may absorb perspiration and other fluids, inspect it for stain and discoloration after each use. When cleaning the cuff, do not machine wash or scrub it. Use a synthetic detergent and gently rub the surface. Air dry thoroughly. Make sure fluid never gets in AIR HOSE.

When storing the unit, do not place heavy objects on it and do not coil AIR HOSE too tightly. When the unit has been stored at a temperature below the freezing point, keep it for at least 1 hour in a warm place before using. Remove the batteries if the instrument is to be stored for an extended period of time. **Keep the batteries out of reach of children.**

We suggest that you have your blood pressure monitor checked every 2 years. This operation may only be performed by the manufacturer or by firms authorized by the manufacturer.

Do not inflate the cuff when it is not wrapped around your arm.

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• WARRANTY

This equipment is guaranteed for the period of 2 years after the date of purchase against manufacturing defects when returned along with the proof of date of purchase to the dealer from whom the purchase was made. During this period, the unit will be repaired or replaced free of charge if the failure is attributable to faulty design or manufacture. This warranty does not cover damage or malfunctions caused by improper handling or use contrary to the instructions in this manual. Please contact your dealer for additional warranty provisions which may remain effective after the manufacturer's warranty period has expired.

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NIHON SEIMITSU SOKKI CO., LTD. 2508-13 Nakago Shibukawa Gunma 377-0293 Japan