NISSEI

MODEL DS-400

DIGITAL BLOOD PRESSURE MONITOR

INSTRUCTIONS ENGLISH

(7-15)
PICTORIALS; 2-6

DIGITALER BLUTDRUCK-MONITOR

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(17-26)
ABBILDUNGEN; 2-6

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CONTENTS

	Page
BASIC PRODUCT INFORMATION	8
PRECAUTIONS FOR USE	9
TIPS ON TAKING YOUR BLOOD PRESSURE	9
BATTERY INSTALLATION	10
CORRECT MEASURING POSTURE	11
WRAPPING THE CUFF	11
MEASUREMENT PROCEDURES	11
MEMORY FUNCTIONS	12
ERROR DISPLAYS AND TROUBLESHOOTING	13
CARE AND MAINTENANCE	15

BASIC PRODUCT INFORMATION

This manual is intended to assist the user in the safe and efficient operation of the Model DS-400 Digital Blood Pressure Monitor. The product must be used in accordance with the procedures contained in this manual and must not be used for purposes other than those described herein. It is important to read and understand the entire manual. In particular, please read carefully and become familiar with the section entitled "TIPS ON TAKING YOUR BLOOD PRESSURE".

Indications for Use

The DS-400 system is intended for the non invasive measurement of systolic and diastolic blood pressure and determination of pulse rate in adult patients, i.e., age 15 and above. The product is recommended for use by patients with labile blood pressure or known hypertension in a home care environment as an adjunct to medical management. The cuff will accommodate an upper arm circumference range of approximately 22 to 32 centimeters (8.7 to 12.6 inches). Pressure is measured over a range of 0 to 300 mmHg and pulse rate over a range of 40 to 160 beats/minute.

The user is cautioned that the system may fail to yield specified measurement accuracy if operated or stored in temperature or humidity conditions outside the limits stated in the specifications section of this manual. Also, an inaccurate reading may result if this instrument is used on a child's arm. Consult your physician if you wish to take a child's blood pressure.

Method of Measurement

This product employs the oscillometric method for measurement of blood pressure and heart rate. The cuff is connected to the main unit and wrapped around the arm. Circuits within the cuff sense the small oscillations in pressure against the cuff produced by the expansion and contraction of the arteries in the arm in response to each heart beat. The amplitude of each pressure waves is measured, converted to millimeters of mercury, and displayed on the LCD as a digital value. A memory circuit stores the 30 most recent measurement results for comparison.

SPECIFICATIONS

Model	: DS-400	Memory	: 30 measurement results, computation of the average of stored data
Operating principle	: Oscillometric method	Operating environment	: +10°C to +40°C, 85% Relative humidity or below
Indicator	: 11 digits LCD display	Storage environment	: -5°C to +50°C, 85% Relative humidity or below
Pressure indicating range	: 0 to 300 mmHg (Cuff pressure)	Cuff	: Coverage arm circumference ; 22 to 32 cm
Measuring range	: 50 to 250 mmHg (systolic) 40 to 180 mmHg (diastolic) 40 to 160 pulse/min. (pulse rate)	Main unit	: Weight ; Approx. 60 g, without battery Size ; 99(W) x 30(D) x 58(H) mm
Accuracy	: 3 mmHg (Cuff pressure) 5% of reading (pulse rate)	Key to symbols	
Inflation	: Squeeze bulb		: Type BF equipment
Deflation	: Constant air release valve		: Important ; Read operating instructions
Exhaust	: Manual rapid exhaust valve		: Segregated disposal required
Power supply	: 1.5 volt "AA" (R6P/LR6) type battery		
Power consumption	: 0.06W (max.)		

Specifications are subject to change without notice due to improvements in performance.

This device complies with EMC (IEC60601-1-2:1993), EN1060-1 and EN1060-3.

PRECAUTIONS FOR USE

Do not use cuffs other than the original cuff contained in this merchandise.

Do not use this instrument without consultation with your doctor if you are under dialysis therapy or on anticoagulants, antiplatelets or steroids. Use of this instrument under such conditions could cause internal bleeding.

To avoid any possibility of accidental strangulation, keep this unit away from children and do not drape AIR HOSE around your neck.

For specific information on your blood pressure, contact your physician. Never make any judgment on measurement results on your own.

Use of this device in areas near portable phones, microwave ovens or other devices with strong electromagnetic field may cause malfunctions.

TIPS ON TAKING YOUR BLOOD PRESSURE

Blood pressure is a measurement of the force exerted by the heart in pumping the blood through the arteries and the resistance by the veins in this flow.

Blood pressure changes easily as influenced by mental and physical factors and is never constant.

In general, blood pressure is highest during the working hours and gradually decreases during the afternoon and evening hours. It is low during sleep and slowly increases at a relatively fast rate after arising from bed.

Causes for Changes in Blood Pressure

- Body movement
- Conversation
- Mental Tension
- Emotions
- Meals
- Drinking Alcohol
- Smoking
- Recent Urination or Bowel Movement
- Temperature
- Changes in the surroundings such as movement or noise, etc.

Blood pressure measured at home tends to be lower than when measured in a hospital, clinic or doctor's office.


When blood pressure is measured in a hospital, it could be higher than when measured at home. This is because you are tense at the hospital and relaxed at home. It is important to know your stable normal blood pressure at home.

Let a qualified physician interpret your blood pressure readings.

Depending on your age, weight and general condition, blood pressure can be slightly different. Consult with your doctor on determining what blood pressure is normal for you.

TIPS ON TAKING YOUR BLOOD PRESSURE

Blood pressure changes shown below.
Upper arm blood pressure date.

180~	(severe)	Hypertension				
160~179	(moderate)					
140~159	(mild)					
130~139	High Normal					
120~129	Normal					
~120	Optimal					
SYS mmHg		80	85	90	100	110
DIA		80	84	89	99	109

Blood pressure classification by WHO (1999)

Rest for approximately five minutes and take your blood pressure while relaxed in a quiet environment.

Take blood pressure measurement with the correct posture and do not move nor speak during measurement.


Avoid exercise, eating, drinking alcohol, smoking and other activities that affect your blood pressure right before a measurement.

Take your blood pressure every day at the same time.

The ambient temperature should be approximately 20°C when you take your blood pressure.

BATTERY INSTALLATION

- 1 Open the battery cover.**
- 2 Install an “AA” type battery into the compartment. Make sure that the polarities correspond to the (+) and (-) marks.**
- 3 Close the battery compartment cover.**

CAUTION! Do not use rechargeable batteries.

battery replacement indication

Replace the battery when the battery replacement indication appears on the display or nothing is displayed even if POWER BUTTON is depressed.

Use alkaline battery to increase battery life. Use of regular Carbon-Zinc battery will necessitate more frequent replacement. *The enclosed battery is for monitoring, and its life may be shorter than that of commercial batteries.*

Since neither the device nor batteries are household waste, follow your national/local recycling rules and dispose of them at appropriate collection sites.

CORRECT MEASURING POSTURE

fig.1

Sit at the table and let the table support your arm as you take the blood pressure measurement. Make sure that the measurement location on the upper arm is at approximately the same height as the heart, and that the forearm is extended naturally on the table and does not move.

You may lie on your back and take the measurement. Look at the ceiling, stay calm, and do not move your neck or body during the measurement. Again, make sure that the measurement location on the arm is at approximately the same height as the heart.

If the cuff is lower (higher) than the heart, the measured reading tends to become larger (smaller).

WRAPPING THE CUFF

fig.2

- 1** With the tube of the cuff positioned in the direction of the fingertips, place it on your left arm.
- 2** Wrap the cuff around the arm with the edge of the cuff approximately 2~3 cm above the elbow.
- 3** Press the surface of the cuff to make sure that it is attached securely.

When wrapping the cuff, wrap it loosely enough around the arm so that two fingers can be placed between the cuff and the arm. If the cuff is wrapped more tightly or loosely than this, inaccurate blood pressure readings may result.

If you are wearing a shirt that might restrict circulation in your upper arm or you roll your sleeve up over the upper arm, the blood flow will be restricted, preventing accurate measurement.

MEASUREMENT PROCEDURES

fig.3

- 1** Insert **TUBE PLUG** into **AIR CONNECTOR**.
- 2** Set **AIR CONTROL VALVE** to **RAPID EXHAUST SIDE**.
- 3** Press **POWER BUTTON**.

"" flashes after the initial display.

Take deep breaths and relax.

Do not move, chat or strain your arm or hand during measurement.

Beep sounds and the inflation mark "▲" starts to flash.

4 Set AIR CONTROL VALVE to CONSTANT AIR RELEASE SIDE.

5 Rapidly squeeze INFLATION BULB to inflate the cuff to 30-40 mmHg above your systolic pressure or if not known to approximately 180 mmHg.

6 Gently place INFLATION BULB on the table.

The displayed value (pressure) decreases.

The pulse is indicated by heart mark "♥" with synchronizing beep.

Long beep sounds and the deflation mark "▼" flashes as the measurement completes.

7 Set AIR CONTROL VALVE to RAPID EXHAUST SIDE and release air from the cuff.

Maximum (SYS:systolic) and minimum (DIA:diastolic) blood pressures and pulse rate (PUL) are displayed.

The measured data will automatically be stored in the memory. *When a measurement resulted in an error, it will not be stored.*

8 Press POWER BUTTON to turn off the power.

If you forget to turn off the unit, it will be automatically turned off after 3 minutes.

Stopping Measurement
Set **AIR CONTROL VALVE** to **RAPID EXHAUST SIDE** and press **POWER BUTTON**.

OVER-PRESSURIZATION

Set **AIR CONTROL VALVE** to **RAPID EXHAUST SIDE** and release air from the cuff.

See page 13.

If the inflation mark flashes again during the measurement, start squeezing **INFLATION BULB** to inflate again to a level about 30 mmHg higher in 4 sec., or the "Err" appears on the display.

See page 13.

Do not take measurements repeatedly for congestion of blood could result in false measurement. Let your arm rest for at least 5 minutes.

MEMORY FUNCTIONS

fig.4

The memory can store up to 30 measurements and their average. When the number of measurements exceeds 30, the oldest data will be replaced by the newest data.

Recalling measured data

1 Press MEMORY BUTTON to see stored data.

When no measurements are stored, nothing will be displayed.

The memory data is displayed for approximately 5 seconds. Approximately 5 seconds after MEMORY BUTTON is released, the apparatus will turn off or return to the current measurement display.

2 Press MEMORY BUTTON to move to next stored data.

Every depression of the button switches among the memory data.

The result with memory data number "1" is the latest stored result and results with larger memory data number indicate the older results.

3 Press POWER BUTTON and turn off the unit.

Storage cannot be cancelled.

The stored data will be retained even after the power is turned off.

The stored data will be lost after battery replacement.

While MEMORY BUTTON is held down for approximately 2 seconds, the memory data numbers will be fast-forwarded to the average of total stored data unless the button is released; the number stops at release of the button.

ERROR DISPLAYS AND TROUBLESHOOTING

ERROR SYMOL / SYMPTOM	CAUSE / CHECK POINT	REMEDY
 MEASURE- MENT ERROR	<p>TUBE PLUG is not correctly inserted.</p> <p>The cuff is not properly wrapped.</p> <p>Measurement could not be made because of moving or talking during measurement.</p> <p>Inflation was not enough.</p>	<p>Make sure that TUBE PLUG is securely inserted.</p> <p>Make sure that the cuff is wrapped properly.</p> <p>Remain still and quiet during measurement.</p> <p>Inflate the cuff when the inflation mark is displayed.</p>
 OVER- PRESSURI- ZATION	<p>The cuff was inflated to the maximum pressure.</p>	<p>Do not inflate the cuff over 300 mmHg.</p>

ERROR DISPLAYS AND TROUBLESHOOTING

ERROR SYMOL / SYMPTOM	CAUSE / CHECK POINT	REMEDY
 IMPROPER OPERATION	POWER BUTTON was accidentally pressed during battery replacement.	Press POWER BUTTON to turn off the power once and press it again to start a measurement.
No display when POWER BUTTON is pressed.	Have the battery polarities been positioned incorrectly? Has the battery run down? Are the battery terminals clean?	Reinsert the battery in the correct position. Replace the battery with new one. Clean the battery terminals with a dry cloth.
No display when MEMORY BUTTON is pressed.	Is there any measured result stored in the memory? Has the battery run down?	Nothing is displayed when there is no measured data stored. The stored data will be lost after battery replacement. Replace the battery with new one.
Cannot make great number of measurements.	Are you using alkaline battery?	Use alkaline battery to increase battery life.
Cannot complete the measurement.	Has the battery run down?	Replace the battery with new one.
Blood pressure is different each time. The reading is extremely low (or high).	Are you measuring with correct posture? Blood pressure readings constantly vary with time of measurement and nervous condition.	Measure with the correct posture. See TIPS ON TAKING YOUR BLOOD PRESSURE.
Pulse rate is too low (or too high).	Did you move during the measurement? Did you make measurement right after exercise?	Remain still during the measurement. Measure again after resting for more than 5 minutes.

ERROR DISPLAYS AND TROUBLESHOOTING

ERROR SYMOL / SYMPTOM	CAUSE / CHECK POINT	REMEDY
The power is automatically turned off.	Have you left the instrument untouched after the measurement?	This does not indicate a problem. The instrument automatically turns off 3 minutes after a measurement or 5 seconds after memory display.

If you cannot get correct measurement with the methods above, do not tamper with the internal mechanism. Contact your dealer. If the unit malfunctions, please return it to the dealer or an authorized service representatives for service according to the warranty.

CARE AND MAINTENANCE

Because the unit includes precision parts, care should be taken to avoid extreme temperature variations, humidity, shock, dust, and direct sunlight. Do not drop or strike the unit. Make sure not to expose the unit to moisture. **This unit is not water resistant.**

Use only a soft, dry cloth to clean the unit. Do not use gasoline, paint thinner, or other strong solvents on the unit. Since the cuff may absorb perspiration and other fluids, inspect it for stain and discoloration after each use. When cleaning the cuff, do not machine wash or scrub it. Use a synthetic detergent and gently rub the surface. Make sure fluid never gets in AIR HOSE. Air dry thoroughly.

When storing the unit, do not place heavy objects on it and do not coil AIR HOSE too tightly. When the unit has been stored at a temperature below the freezing point, keep it for at least 1 hour in a warm place before using it. Remove the battery if the instrument is to be stored for an extended period of time. **Keep the batteries out of reach of children.**

We suggest that you have your blood pressure monitor checked every 2 years. This operation may only be performed by the manufacturer or by firms authorized by the manufacturer.

Do not inflate the cuff when it is not wrapped around your arm.

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WARRANTY

This equipment is guaranteed for the period of 2 years after the date of purchase against manufacturing defects when returned along with the proof of date of purchase to the dealer from whom the purchase was made. During this period, the unit will be repaired or replaced free of charge if the failure is attributable to faulty design or manufacture. This warranty does not cover damage or malfunctions caused by improper handling or use contrary to the instructions in this manual. Please contact your dealer for additional warranty provisions which may remain effective after the manufacturer's warranty period has expired.

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Manufacturer: **NIHON SEIMITSU SOKKI CO., LTD.**

2508-13 Komochi Kitagunma Gunma Japan

EC-Representative: **JPI Trading GmbH**

Holzbrücke 7, 20459 Hamburg, Germany ☎ 040-82 22 80 510

web site <http://www.nissei-jp.com>

A113343-1-A