WRIST TYPE DIGITAL BLOOD PRESSURE MONITOR WS-1300

NISSEI CE 0123

Manufacturer:

NIHON SEIMITSU SOKKI CO., LTD. 2508-13 Nakago Shibukawa Gunma 377-0293 Japan

EC-Representative:

Nissei Healthcare (UK) Ltd. Henfield, BN5 9SJ UK

web site <http://www.nisseihealthcare.com>

INSTRUCTIONS

ENGLISH

→ page 9

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DEUTSCH

→ Seite 20

INSTRUCTIONS

FRENCH

→ page 32

ISTRUZIONI

ITALIANO

→ pagina 44

INSTRUCCIONES

ESPAÑOL

→ página 56

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JĘZYK POLSKI

→ str 68

PICTORIALS

ABBILDUNGEN

IMAGES

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FIGURAS

ILUSTRACJE

Fig.1

PARTS NAMES AND PRODUCT COMPONENTS

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

8

9

- ① DISPLAY
- ② POWER BUTTON
- ③ MEMORY M1 BUTTON
- ④ MEMORY M2 BUTTON
- ⑤ BATTERY COVER
- ⑥ WRIST CUFF
- ⑦ CARRYING CASE
- ⑧ 2 LR03 (AAA) ALKALINE BATTERIES
- ⑨ INSTRUCTION MANUAL

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Fig.2

**INSTALLING BATTERIES AND
SETTING THE CLOCK**

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str. 71










Fig.3

**WRAPPING THE CUFF
AROUND YOUR WRIST**

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
BLOOD PRESSURE
CLASSIFICATION BY WHO
(1999)

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
CONTENTS

PICTORIALS	2	CORRECT MEASURING POSTURE	14
BASIC PRODUCT INFORMATION	9	MEASURING METHOD	14
PRECAUTIONS FOR USE	11	MEMORY FUNCTION	16
TIPS ON TAKING YOUR BLOOD PRESSURE	11	ERROR DISPLAYS AND TROUBLESHOOTING	17
INSTALLING BATTERIES AND SETTING THE CLOCK	12	CARE AND MAINTENANCE	18
WRAPPING THE CUFF AROUND YOUR WRIST	13	TECHNICAL DESCRIPTION	19

BASIC PRODUCT INFORMATION

This manual is intended to assist the user in the safe and efficient operation of the Model WS-1300 Wrist Type Digital Blood Pressure Monitor. The product must be used in accordance with the procedures contained in this manual and must not be used for purposes other than those described herein. It is important to read and understand the entire manual. In particular, please read carefully and become familiar with the section entitled “TIPS ON TAKING YOUR BLOOD PRESSURE”.

The user is cautioned that the system may fail to yield specified measurement accuracy if operated or stored in temperature or humidity conditions outside the limits stated in the specifications

section of this manual.

INDICATIONS FOR USE

The WS-1300 system is intended for the non invasive measurement of systolic and diastolic blood pressure and determination of pulse rate in adult, i.e., age 15 and above. The product is recommended for use by people with labile (unstable) blood pressure or known hypertension in a home care environment as an adjunct to medical management. The cuff will accommodate a wrist circumference range of approximately 125 to 215 millimeters. Pressure is measured over a range of 0 to 300 mmHg and pulse rate over a range of 40 to 160 beats/minute.

METHOD OF MEASUREMENT

This product employs the oscillometric method for measurement of blood pressure and pulse rate. The entire system, including the cuff, is a single wrist-mounted unit. The cuff is wrapped around the wrist and automatically pressurized. Circuits within the cuff sense the small oscillations in pressure against the cuff produced by the expansion and contraction of the arteries in the wrist in response to each heart beat. The amplitude of each pressure waves is measured, converted to millimeters of mercury, and displayed on the LCD as a digital value. Responding to heart beats, the electromagnetic deflation control valve maintains the constant deflation rate irrespective of differences in wrist size. A memory circuit stores the 60 most recent measurements for comparison.

SPECIFICATIONS

Model	: WS-1300
Operating Principle	: Oscillometric method
Indicator	: 12 digit Liquid Crystal Display
Pressure Indication Range	: 0 to 300 mmHg (cuff pressure)
Measuring Range	: 50 to 250 mmHg (systolic) : 40 to 180 mmHg (diastolic) : 40 to 160 pulses/min (pulse rate)
Accuracy	: ±3 mmHg (cuff pressure) : ±5 % of reading (pulse rate)

Inflation	: Automatic (air pump)
Deflation	: Automatic (electric control valve)
Exhaust	: Automatic quick exhaust valve
Power Supply	: 2 AAA alkaline (LR03) dry cells
Memory	: 60 measurements x 2 and averages
Applied Pressure	: 190 mmHg (fixed)
Operating Environment	: +10°C to+40°C, 85% relative humidity or below
Storage Environment	: -5°C to+50°C, 85% relative humidity or below
Cuff	: Pre-formed type cuff
Wrist circumference	: Approx. 125 to 215 mm
Main Unit	: Size ; 70(W) x 70(H) x 29.9(D) mm : Weight ; Approx. 110 g., without batteries
Key to symbols	:  Type BF equipment :  Important ; Read operating instructions
Classification	: Internally powered equipment

Specifications are subject to change without prior notice due to improvements in performance and quality.
This device complies with EMC (IEC60601-1-2, EN60601-1), EN1060-1 and EN1060-3.

PRECAUTIONS FOR USE

Do not use this instrument without consultation with your doctor if you are under dialysis therapy or on anticoagulants, antiplatelets or steroids. Use of this instrument under such conditions could cause internal bleeding.

For specific information on your blood pressure, contact your physician. Never make any judgment on your own

regarding measurement results.

Use of this device in areas near portable phones, microwave ovens or other devices with strong electromagnetic field may cause malfunctions.

TIPS ON TAKING YOUR BLOOD PRESSURE

Blood pressure is the force exerted by the heart in pumping the blood through the arteries and the resistance by the veins to this flow.

Blood pressure varies all the time, influenced by mental and physical factors and is never constant.

In general, blood pressure is highest during the working hours and gradually decreases during the afternoon and evening hours. It is low during sleep and increases at a relatively fast rate after arising from bed.

Causes for Changes in Blood Pressure

- Body movement
- Emotions
- Recent Urination or Bowel Movement
- Changes in the surroundings such as movement or noise, etc.
- Conversation
- Eating
- Drinking Alcohol
- Mental Tension
- Smoking
- Temperature

Blood pressure measured at home tends to be lower than when measured in a hospital, clinic or doctor's office.

This is because you are tense at the hospital and relaxed at home. It is important to know your stable normal blood pressure at home.

Blood pressure at the wrist may be different from that at the upper arm.

As the thickness of the blood vessels differs depending on measuring position, blood pressure value measured at the wrist may be slightly different from that measured at the upper arm. In most cases, the difference for a healthy person will be within about ± 10 mmHg for the systolic and the diastolic blood pressure. Caution is required for the following persons: persons with hypertension, diabetes, liver trouble, hardened arteries, peripheral hindrance of

circulation, etc., as there may be a difference between the values measured at the wrist and those measured at the upper arm.

Let a qualified physician interpret your blood pressure readings.

Depending on your age, weight and general condition, blood pressure can be slightly different. Consult with your doctor on determining what blood pressure is normal for you.

Before taking blood pressure, rest for approximately five minutes and take your blood pressure while relaxed in a quiet environment.

Measure blood pressure using the correct posture and do not move nor speak during measurement.

Avoid exercise, eating, drinking alcohol, smoking and other activities that affect your blood pressure right before a measurement.

Take your blood pressure at the same time every day.

The ambient temperature should be approximately 20°C when you take your blood pressure.

INSTALLING BATTERIES AND SETTING THE CLOCK

Fig.2 *page 4*

1. Remove the battery cover in the arrow direction, lightly pressing its sides.
2. Install two AAA alkaline dry cells in the battery compartment, giving attention to their polarity.
The cells can be easily installed by pushing them while pressing their (-) end against the spring.
3. Replace the battery cover.
4. Year flashes on the display.

Clock can be set only after the batteries are reinstalled. Setting the clock correctly ensures storage of the measurement with the correct date and time.

Use memory M1 button to increase the flashing year and memory M2 button to decrease the year. Press POWER BUTTON to confirm and move to next step.

5. Month flashes. Adjust with M1 and M2 buttons and press POWER

Do not force the battery cover into position.

BUTTON to confirm.

6. Day flashes. Adjust with M1 and M2 buttons and press POWER BUTTON to confirm.

7. Hour flashes. Adjust with M1 and M2 buttons and press POWER BUTTON to confirm.

The clock operates in 24-hour mode.

8. Minute flashes. Adjust with M1 and M2 buttons and press POWER BUTTON to confirm.

The clock is set and the unit is turned off.

BATTERY REPLACEMENT INDICATION

Replace all the batteries when the battery replacement indication

appears on the display or nothing is displayed when POWER BUTTON is depressed.

About 300 measurements can be made after installation of new AAA alkaline dry cells. *The enclosed batteries are for monitoring, and their life may be shorter than that of commercial batteries.*

 The used electrical and electronic products are not household waste. Follow your national/local recycling rules to dispose of them properly. In the EU countries, please refer to waste management symbol(s) marked on the package or the instrument.

Do not use rechargeable batteries.

WRAPPING THE CUFF AROUND YOUR WRIST

Fig.3 page 5

1. With the palm of your left hand up, put the cuff on your wrist so that the display screen is on the same side as your palm. If it is impossible to put the cuff on your left wrist, use your right wrist for measurement.

2. Adjust the cuff, with your palm up, until its top edge is positioned 5 to 10 mm from the bottom of your palm. Place the unit at the

center of your wrist.

3. Fasten the cuff around your wrist so that there is no space between the cuff and your wrist. The cuff should fit snugly. Press the surface of the cuff to make sure the hook and loop attaches securely. Position the cuff directly against the skin. Take care that your clothing is not caught beneath the cuff.

CORRECT MEASURING POSTURE

Fig.4 page 5

1. Sit on a chair.
2. Lightly open your left hand, with your palm up, and put your elbow on a table or stand.
3. Adjust the cuff to the same height as your heart with a comfortable support such as a folded towel etc.

WHEN THERE IS NO TABLE AVAILABLE

Sit on chair. Adjust the cuff to the same height as your heart by lightly applying your left arm to your breast. Lightly support your arm by your right hand during measurement, but do not touch the cuff.

TAKING A MEASUREMENT IN BED

Lie on your back. Adjust the cuff to the same height as your heart with a comfortable support such as a folded towel etc.

Measured data may vary slightly depending on the posture during measurement.

Measurements should always be taken using the same wrist in the same position.

If the cuff is lower (higher) than the heart, the measured reading tends to become larger (smaller).

MEASURING METHOD

Fig.5 page 6

Breathe deeply and relax!

Make sure the cuff is wrapped correctly and it is at the same height as your heart.

Do not move, chat or strain your arm or hand during measurement.

1. Press POWER BUTTON.

The cuff is automatically inflated to the applied pressure.

TO STOP MEASUREMENT

Press POWER BUTTON and the unit stops inflation, discharges air rapidly, and then turns off.

Pressure (displayed value) starts to decrease and pulse is indicated by heart mark.

When the measurement is complete, air is automatically released from the cuff.

Upper and lower blood pressures and pulse rate are displayed.

Heart mark flashes when pulse rhythm detected during measurement was irregular.

2. Press either one of memory buttons, M1 or M2, and the result is stored in the selected bank.

The selected memory bank number is shown on the display. The result is stored in the bank selected when the unit is turned off. When a measurement resulted in an error, it will not be stored.

3. Press POWER BUTTON to turn off the power.

If you forget to turn off the unit, it will automatically turn off after 3 minutes.

Do not execute repeated measurements. This will cause engorgement of the wrist and correct values will not be obtained. Let your wrist rest for at least 5 minutes.

IRREGULAR PULSE RHYTHM INDICATION

Heart mark flashing in measurement result display indicates irregular pulse rhythm.

Irregular pulse rhythm can be a result of body movement during measurement or can be because of arrhythmias. Although frequent appearance of the indication under quiet measurements may suggest arrhythmias, do not make any judgment on your own before consulting with your doctor.

An example of irregular pulse rhythm compared to regular pulse rhythm

AUTOMATIC REPRESSURIZATION

If the pressure applied is judged insufficient in the early phase of measurement or if movement of the hand or wrist occurs during measurement, the unit will inflate again to a level about 40 mmHg higher. Automatic Repressurization is repeated until a measurement is made. However, this does not indicate a problem. To manually control pressurization, hold down POWER BUTTON immediately after the start of inflation and release it when the pressure reaches the level 50 mmHg higher than the expected maximum systolic pressure. If the pressure value has exceeded 190 mmHg, the inflation will stop when the key is released. Pressure can be increased to approximately 300 mmHg.

Measured result is stored in either M1 or M2. Each of two memory banks can hold up to 60 results and their average. When the number of measurements exceeds 60, the oldest data will be deleted to record new data. You may decide which bank to store your measurement results to avoid data mixture with someone else's, or may use two banks to save data measured in the morning and evening separately.

RECALLING STORED DATA

1. Press memory M1 button or memory M2 button to see stored data.

The average of the stored result in the selected bank is displayed.

The latest result stored is displayed when there is only a single result and the clock display remains when there are no results stored.

The memory data is displayed for approximately 30 seconds. Approximately 30 seconds after memory button is released, the apparatus will turn off.

2. Press memory button to move to next data.

The memory data number 1 is the latest among the stored data in the selected memory bank.

The indication at the top of the display alternately changes from memory data number to date and to time.

Heart mark flashes when pulse rhythm detected during measurement was irregular.

3. Every depression of memory button switches among the memory data.

Readings are recalled from the latest reading to the oldest.

DELETING STORED DATA

1. Show the individual result to be erased or the average of a memory bank to clear all the memory bank data.

2. Press and hold down either one of memory buttons: M1 button to erase specified data in memory bank M1, and M2 button to erase specified data in memory bank M2.

The displayed data starts flashing.

3. Hold down the button until nothing but the memory bank number is displayed.

ERROR DISPLAYS AND TROUBLESHOOTING

OVER-PRESSURIZATION

The pressurization was not enough to make a measurement even though the cuff was inflated to the maximum value because of wrist movement or talking during measurement. **Do not move during measurement.**

MEASUREMENT ERROR

Measurement could not be made because of moving or talking during measurement. **Remain still and quiet during measurement.**

INFLATION ERROR

The cuff is not properly applied. **Properly apply the cuff.**

If Err-2 appears after re-trying to apply the cuff properly, assume the cuff is damaged. **Contact your dealer.**

DEFLATION ERROR

Movement or talking was detected during measurement. **Turn off the power and make another measurement while being quiet and still.**

IMPROPER OPERATION

POWER BUTTON was accidentally pressed during battery replacement. **Press POWER BUTTON once to turn off the power and press it again to start a measurement.**

No display when POWER BUTTON is pressed.

Batteries are not correctly inserted. **Reinsert the batteries in the correct position.**

Battery terminals are not clean. **Clean with a dry cloth.**

Batteries are exhausted. **Replace all batteries with new ones.**

The reading is extremely low (or high). Blood pressure is different each time.

Blood pressure readings constantly vary with time of measurement and both mental and physical conditions. **See TIPS ON TAKING YOUR BLOOD PRESSURE.**

Measurement is not taken with correct posture or the cuff is wrapped incorrectly. **Reconfirm the measurement procedures.**

The reading was effected from moving or talking. **Remain**

still and quiet during measurement.

Measurement was taken right after exercise or other activities influencing the reading. **Measure again after resting for more than 5 minutes. See TIPS ON TAKING YOUR BLOOD PRESSURE.**

If you cannot get correct measurement with the methods above, stop using the unit and contact your dealer. Do not tamper with the internal mechanism. If the unit malfunctions, please return it to the dealer or an authorized service representative for service according to the warranty.

CARE AND MAINTENANCES

Because the unit includes precision parts, care should be taken to avoid extreme temperature variations, humidity, shock, dust, and direct sunlight. Do not drop or strike the unit. Make sure not to expose the unit to moisture. **This unit is not water resistant.**

Clean the body with a cloth moistened with water and a mild detergent, and then wipe dry. Do not use gasoline, paint thinner, or other strong solvents. Since the cuff may absorb perspiration and other fluids, inspect it for stain and discoloration after each use. When cleaning the cuff, use a synthetic detergent and gently rub the surface.

When storing the unit, do not place heavy objects on it. When the unit has been stored at a temperature below the freezing point, keep it for at least 1 hour in a warm place before using it.

Remove the batteries if the instrument is to be stored for an extended period of time. **Keep the batteries out of reach of children.**

We suggest that you have your monitor checked every 2 years. This operation may only be performed by the manufacturer or by firms authorized by the manufacturer.

Do not inflate the cuff (start a measurement) when it has not been wrapped around your wrist. Do not try to stretch or twist the cuff, and do not bring it into contact with sharp or pointed objects.

Do not disassemble or modify the unit.

TECHNICAL DESCRIPTION

WS-1300 complies with the EMC, electromagnetic compatibility, standard, IEC60601-1-2. Refer to the tables below for specific information regarding compliance to the standard.
 WS-1300, as a medical electrical equipment, needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided below.
 Portable and mobile RF communications equipments can affect the device.
 The use of accessories other than those specified in this manual may result in increased emissions or decreased immunity of the device.
 WS-1300 should not be used adjacent to or stacked with other equipment.

Table 201 - Guidance and manufacturer's declaration - electromagnetic emissions -

WS-1300 is intended for use in the electromagnetic environment specified below. The customer or the user of WS-1300 should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	WS-1300 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	WS-1300 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	N/A	

Table 204 - Guidance and manufacturer's declaration - electromagnetic immunity -

WS-1300 is intended for use in the electromagnetic environment specified below. The customer or the user of WS-1300 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2,5 GHz	N/A 3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of WS-1300, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance N/A $d = 1.2\sqrt{P}$, 80 MHz to 800 MHz $d = 1.2\sqrt{P}$, 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

^a Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radios broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which WS-1300 is used exceeds the applicable RF compliance level above, WS-1300 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating WS-1300.
^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V] V/m.

Table 202 - Guidance and manufacturer's declaration - electromagnetic immunity -

WS-1300 is intended for use in the electromagnetic environment specified below. The customer or the user of WS-1300 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	N/A	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	N/A	N/A
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	N/A	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_i (>95% dip in U_i) for 0,5 cycle 40% U_i (60% dip in U_i) for 5 cycles 70% U_i (30% dip in U_i) for 25 cycles <5% U_i (>95% dip in U_i) for 5 sec	N/A	N/A
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_i is the a.c. mains voltage prior to application of the test level.

Table 206 - Recommended separation distances between portable and mobile RF communications equipment and WS-1300 -

WS-1300 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of WS-1300 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and WS-1300 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter, W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz, N/A	80 MHz to 800 MHz, $d = 1.2\sqrt{P}$	800 MHz to 2,5 GHz, $d = 1.2\sqrt{P}$
0,01	N/A	0,12	0,23
0,1	N/A	0,38	0,73
1	N/A	1,2	2,3
10	N/A	3,8	7,3
100	N/A	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
 NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
 NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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